



"ISO 9001:2008 & EN ISO 13485:2012"
CERTIFIED COMPANY
"CE MARKED PRODUCTS"

DRY BATH INCUBATOR

WITH INTERCHANGEABLE BLOCK



Model No : DBI - 33

NOTE :- The Heating Block can be replaced by user as they need.

SPECIFICATION :-


Mains Input : 230V/50Hz | 110V/60Hz A.C. Supply
W x D x H : 168 (W) x 232 (D) x 93 (H)
Body : Metal Body
Weight : 2.3 Kg. Net
Temp. Controller : Digital Type

Temp. Range : 37°C To 70°C
Accuracy : $\pm 0.5^\circ\text{C}$ @ 37°C, $\pm 3\%$ Except 37°C

HEATING BLOCK : INTERCHANGEABLE

Timer : 1 To 999 (Digital Type)

Block - 1



12.7 MM 24 HOLE

Block - 2



10.5 MM 28 HOLE

Block - 3



**18.5 MM - 4 HOLE
12.7 MM - 6 HOLE
15.5 MM - 5 HOLE
10.5 MM - 7 HOLE**